



## Tulip Power Injector (PI4) - Injection Volume Chart

<b>Small Barrel</b>	<b>Per gun click</b>	<b>Per full pull (4 clicks)</b>
10cc BD Luer	0.2ml	0.8ml
10/12cc Monoject Luer Lock	0.25ml	1.0ml
10cc Monoject for Oral Medication	0.25ml	1.0ml
20cc BD Luer Lock	0.4ml	1.6ml
20cc Monoject Luer Lock	0.42ml	1.7ml

  

<b>Large Barrel</b>	<b>Per gun click</b>	<b>Per full pull (4 clicks)</b>
60cc Toomey Tip	0.8ml	3.2ml
60cc Nipro Luer Lock	1.1ml	4.4ml
60cc Monoject Slip Tip	0.92ml	3.7ml